



Product Datasheet

ACK1 (phospho-Y284) antibody (orb304637)


Description	Rabbit polyclonal antibody to ACK1 (phospho-Y284)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding Y284 of human ACK1 protein. The exact sequence is proprietary.
Target	TNK2
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	O54967 , Q07912 , Q5U2X5
Entrez	10188 , 51789 , 303882
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of ACK1 (Phospho-Y284)



Immunohistochemical analysis of ACK1 (Phospho-Y284)



Immunofluorescent analysis of ACK1 (Phospho-Y284)