



Product Datasheet

ZIC1/2/3 antibody (orb304668)


Description	Rabbit polyclonal antibody to ZIC1/2/3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human ZIC1/2/3. The exact sequence is proprietary.
Target	ZIC1; ZIC2; ZIC3
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q62520 , O60481 , Q15915 , P46684 , Q62521 , O95409
Entrez	22773 , 7545 , 7546 , 22771 , 7547
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of ZIC1/2/3 expres...



Immunohistochemical analysis of ZIC1/2/3...



Immunofluorescent analysis of ZIC1/2/3 S...