



Product Datasheet

PRAS40 (phospho-T246) antibody (orb304682)

Description	Rabbit polyclonal antibody to PRAS40 (phospho-T246)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T246 of human PRAS40 protein. The exact sequence is proprietary.
Target	AKT1S1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q9D1F4, Q96B36
Entrez	67605, 84335
Dilution Range	WB: 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of PRAS40 (Phospho-T246)



Immunohistochemical analysis of PRAS40 (Phospho-T246)