



Product Datasheet

Ataxin 1 (phospho-S775) antibody (orb304691)


Description	Rabbit polyclonal antibody to Ataxin 1 (phospho-S775)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S775 of human Ataxin 1 protein. The exact sequence is proprietary.
Target	ATXN1
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q63540 , P54254 , P54253
Entrez	25049 , 6310 , 20238
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Ataxin 1 (Phospho-S775)



Immunohistochemical analysis of Ataxin 1 (Phospho-S775)



Immunofluorescent analysis of Ataxin 1 (Phospho-S775) (...