

Product Datasheet

CDC6 (phospho-S54) antibody (orb304718)

Description Rabbit polyclonal antibody to CDC6 (phospho-S54)

Species/Host Rabbit

Reactivity Human, Mouse, Primate, Rat

Conjugation Unconjugated

Tested Applications IH, WB

Immunogen KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S54 of human CDC6 protein. The exact sequence is proprietary.

Target CDC6

Preservatives Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note For research use only

Clonality Polyclonal


Source Rabbit

Uniprot ID [Q99741](#), [O89033](#)


Entrez [23834](#), [990](#)

Dilution Range WB: 1:500-1000, IHC-P: 1:100-200

Expiration Date 12 months from date of receipt.



Western blot analysis of CDC6 (Phospho-S54)



Immunohistochemical analysis of CDC6 (Phospho-S54)