



Product Datasheet

C/EBP beta (phospho-T235) antibody (orb304722)


Description	Rabbit polyclonal antibody to C/EBP beta (phospho-T235)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IF, IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding T235 of human C/EBP beta protein. The exact sequence is proprietary.
Target	CEBPB
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	P17676, P28033, P21272
Entrez	24253, 1051, 12608
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200, IF/ICC: 1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of C/EBP beta (Pho...)



Immunohistochemical analysis of C/EBP be...



Immunofluorescent analysis of C/EBP beta...