


Product Datasheet

IKK alpha/beta (phospho-S176/177) antibody (orb304723)

Description

Rabbit polyclonal antibody to IKK alpha/beta

kDa A B C



Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Primate, Rat, Zebrafish

Conjugation

Unconjugated


Tested Applications

IH, WB

Western blot analysis of IKK alpha/beta ...

Immunogen

KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S176/177 of human IKK alpha/beta protein. The exact sequence is proprietary.



Immunohistochemical analysis of IKK alpha/beta ...

Target

CHUK; IKBKB

Preservatives

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Form/Appearance

Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.

Storage

Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Note

For research use only

Clonality

Polyclonal

Source

Rabbit

Uniprot ID

[O14920](#), [O88351](#), [Q9QY78](#), [O15111](#), [Q60680](#)

Entrez

[1147](#), [16150](#), [84351](#), [3551](#)

Dilution Range

WB: 1:500-1000

Expiration Date

12 months from date of receipt.