



Product Datasheet

DAXX (phospho-S668) antibody (orb304735)


Description	Rabbit polyclonal antibody to DAXX (phospho-S668)
Species/Host	Rabbit
Reactivity	Human, Mouse, Primate, Rat
Conjugation	Unconjugated
Tested Applications	IH, WB
Immunogen	KLH-conjugated synthetic phosphopeptide corresponding to residues surrounding S668 of human DAXX protein. The exact sequence is proprietary.
Target	DAXX
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q8VIB2 , O35613 , Q9UER7
Entrez	13163 , 1616
Dilution Range	WB: 1:500-1000, IHC-P: 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of DAXX (Phospho-S668) antibody.



Immunohistochemical analysis of DAXX (Phospho-S668) antibody.



Direct ELISA antibody dose-response curve for DAXX (Phospho-S668) antibody.