

Product Datasheet

DNA Polymerase gamma 2 peptide (orb304752)

Description	The peptide is used to DNA Polymerase gamma 2 antibody ractivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	POLG2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9QZM2, Q9UHN1
Entrez	50776, 11232
Expiration Date	12 months from date of receipt.