

Product Datasheet

R3HCC1L peptide (orb304777)

Description	The peptide is used to R3HCC1L antibody ractivity
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	R3HCC1L
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q7Z5L2, Q8BJM3
Entrez	52013, 27291
Expiration Date	12 months from date of receipt.