

Product Datasheet

RAB11FIP1 peptide (orb304778)

Description	The peptide is used to RAB11FIP1 antibody ractivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	RAB11FIP1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q3B7T9, Q6WKZ4, Q9D620
Entrez	80223, 75767, 498650
Expiration Date	12 months from date of receipt.