

Product Datasheet

NF-kappaB p65 (phospho-T435) peptide (orb304789)

Description	The peptide is used to NF-kappaB p65 (phospho-T435) antibody ractivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	RELA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q04207 , Q04206
Entrez	19697 , 5970
Expiration Date	12 months from date of receipt.