

Product Datasheet

alpha synuclein (phospho-S129) peptide (orb304841)

Description	alpha synuclein (phospho-S129) peptide
Reactivity	Human, Mouse, Primate, Rat, Sheep
Tested Applications	Blocking
Target	SNCA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P37840, P37377, O55042
Entrez	29219, 6622, 20617
Expiration Date	12 months from date of receipt.