

Product Datasheet

TCEAL8 peptide (orb304865)

Description	The peptide is used to TCEAL8 antibody ractivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	TCEAL8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8IYN2, Q9CZY2, Q6I7R5
Entrez	90843, 66684, 367909
Expiration Date	12 months from date of receipt.