

Product Datasheet

PIAS2 antibody (orb30489)

Description	Rabbit polyclonal antibody to PIAS2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	PIAS2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	68240
Uniprot ID	075928
NCBI	NP_004662.2, NP_775298.1
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.

