

Product Datasheet

CaMK2 beta/gamma (phospho-T287) peptide (orb304948)

Description	The peptide is used to CaMK2 beta/gamma (phospho-T287) antibody ractivity
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	CAMK2B; CAMK2G
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P28652, P08413, P11730, Q13555, Q13554, Q923T9
Entrez	816, 12323, 12325, 818, 24245, 171140
Expiration Date	12 months from date of receipt.