

Product Datasheet

c-CBL (phospho-Y700) peptide (orb304951)

Description	The peptide is used to c-CBL (phospho-Y700) antibody ractivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CBL
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P22681, P22682
Entrez	12402, 867
Expiration Date	12 months from date of receipt.