

Product Datasheet

IKK alpha/beta (phospho-S176/177) peptide (orb304964)

Description	The peptide is used to IKK alpha/beta (phospho-S176/177) antibody reactivity
Reactivity	Human, Mouse, Porcine, Primate, Rat, Zebrafish
Tested Applications	Blocking
Target	CHUK; IKBKB
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q60680 , O88351 , O15111 , Q9QY78 , O14920
Entrez	1147 , 16150 , 3551 , 84351
Expiration Date	12 months from date of receipt.