


Product Datasheet

CD193 antibody (orb305468)

Description	Mouse monoclonal antibody to CD193
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	human CD193 transfectants
Target	CD193
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P51677
Entrez	1232
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...