



Product Datasheet

CD154 antibody (orb305469)

Description	Mouse monoclonal antibody to CD154
Reactivity	Human, Primate
Conjugation	Unconjugated
Tested Applications	FC, ICC
Immunogen	human CD154 fusion protein
Target	CD154
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-6 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P27548
Entrez	21947
Dilution Range	Flow cytometry: Recommended dilution: 2-6 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD154 positive cells...



Flow cytometry surface staining pattern ...