

Product Datasheet

DOT1L antibody (orb305479)

Description	Rabbit polyclonal antibody to DOT1L
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	DOT1L
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	164856
Uniprot ID	Q8TEK3
NCBI	NP_115871.1
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.

