

Product Datasheet

BNIP3L antibody (orb305691)

Description	Rabbit polyclonal antibody to BNIP3L
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	BNIP3L
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23930
Uniprot ID	O60238
NCBI	NP_004322.1
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.

