

Product Datasheet

PDGFR beta antibody (orb305712)

Description	Rabbit polyclonal antibody against PDGFR beta
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	PDGFRB
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	123968
Uniprot ID	P09619
NCBI	NP_002600.1
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.

