



Product Datasheet

BRD4 antibody (orb305724)

Description	Rabbit polyclonal antibody to BRD4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic Peptide
Target	BRD4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	152219
Uniprot ID	O60885
NCBI	NP_055114.1 , NP_490597.1
Dilution Range	IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ca...