

Product Datasheet

PTK2 antibody (orb305742)

Description	Mouse monoclonal antibody to PTK2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	PTK2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	119233
Uniprot ID	Q05397
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

