




Product Datasheet

GATA3 antibody (orb305757)

Description	Mouse monoclonal antibody to GATA3
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Recombinant Protein
Target	GATA3
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2b, κ
Clonality	Monoclonal
MW	47916
Uniprot ID	P23771
Dilution Range	IF/ICC: 1:25, FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cel...



Western blot analysis of CEM, HL-60, SH-...