



Product Datasheet

IL4R antibody (orb305783)

Description	Rabbit polyclonal antibody to IL4R
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	IL4R
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	89658
Uniprot ID	P24394
Dilution Range	FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T, A431 cell...