




Product Datasheet

AK4 antibody (orb305785)

Description	Rabbit polyclonal antibody to AK4
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	AK4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25268
Uniprot ID	P27144
Dilution Range	IF/ICC: 1:25, IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



IHC-P staining of mouse kidney section t...



Immunohistochemical staining of human ki...