



Product Datasheet

MAFK antibody (orb305786)

Description	Mouse monoclonal antibody to MAFK
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Recombinant Protein
Target	MAFK
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	17523
Uniprot ID	O60675
Dilution Range	FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lin...