

Product Datasheet

MAGED1 antibody (orb305789)

Description	Mouse monoclonal antibody to MAGED1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	MAGED1
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	86161
Uniprot ID	Q9Y5V3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

