



Product Datasheet

PPARA antibody (orb305792)

Description	Mouse monoclonal antibody to PPARA
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Recombinant Protein
Target	PPARA
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	52225
Uniprot ID	Q07869
Dilution Range	IHC-P: 1:25, IF/ICC: 1:25, FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Immunohistochemical staining of human sk...