




Product Datasheet

FBXL17 antibody (orb305799)

Description	Rabbit polyclonal antibody to FBXL17
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	FBXL17
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	75695
Uniprot ID	Q9UF56
Dilution Range	IHC-P: 1:25, FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



IHC-P staining of human kidney section t...



Immunohistochemical staining of rat kidn...