




Product Datasheet

ADRA1D antibody (orb305807)

Description	Rabbit polyclonal antibody to ADRA1D
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ADRA1D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60463
Uniprot ID	P25100
Dilution Range	IHC-P: 1:25, FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



IHC-P staining of human kidney section t...



Immunohistochemical staining of human sm...