




Product Datasheet

FAM216B antibody (orb305830)

Description	Rabbit polyclonal antibody to FAM216B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	FAM216B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	16362
Uniprot ID	Q8N7L0
Dilution Range	IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of MCF-7, SK-BR-3,...