



Product Datasheet

AURKA antibody (orb305839)

Description	Mouse monoclonal antibody to AURKA
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	AURKA
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,k
Clonality	Monoclonal
MW	45823
Uniprot ID	O14965
Dilution Range	IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HT-29 cell line...