




Product Datasheet

vWF antibody (orb305844)

Description	Mouse monoclonal antibody to vWF
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Recombinant Protein
Target	VWF
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	309265
Uniprot ID	P04275
Dilution Range	FACS: 1:100, IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of VWF antibody