




Product Datasheet

SHB antibody (orb305848)

Description	Rabbit polyclonal antibody to SHB
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	SHB
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55042
Uniprot ID	Q15464
Dilution Range	IHC-P: 1:25, FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of HepG2, Jurkat, ...