




Product Datasheet

GLI2 antibody (orb305854)

Description	Rabbit polyclonal antibody to GLI2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Target	GLI2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	167783
Uniprot ID	P10070
Dilution Range	IF/ICC: 1:25, FACS: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HeLa cell...



Western blot analysis of 293, Jurkat cel...