




Product Datasheet

CAPN2 antibody (orb305855)

Description	Mouse monoclonal antibody to CAPN2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Target	CAPN2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	79995
Uniprot ID	P17655
Dilution Range	IHC-P: 1:25, FACS: 1:100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot analysis of A549, A431, SK-...