

Product Datasheet

HSP90AA1 antibody (orb305874)

Description	Rabbit polyclonal antibody to HSP90AA1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	HSP90AA1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	84660
Uniprot ID	P07900
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

