




Product Datasheet

SMAD1 antibody (orb305878)

Description	Rabbit polyclonal antibody to SMAD1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Target	SMAD1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52260
Uniprot ID	Q15797
Dilution Range	IF/ICC: 1:25, IHC-P: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human he...



Immunohistochemical staining of paraffin...