



Product Datasheet


TUFM antibody (orb305889)

Description	Rabbit polyclonal antibody to TUFM
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	TUFM
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49542
Uniprot ID	P49411
Dilution Range	IF/ICC: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 cell line

...



Western blot analysis of HepG2, Hela, mo...