



Product Datasheet

SNRPG antibody (orb305987)

Description	Rabbit polyclonal antibody to SNRPG
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Target	SNRPG
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	8496
Uniprot ID	P62308
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Human Cer...