

Product Datasheet

PIK3R1 antibody (orb306009)

Description	Rabbit polyclonal antibody to PIK3R1
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	PIK3R1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	83598
Uniprot ID	P27986
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

