

Product Datasheet

SYCP3 antibody (orb306017)

Description	Rabbit polyclonal antibody to SYCP3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	SYCP3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	27729
Uniprot ID	Q8IZU3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

