

Product Datasheet

HIST1H3A antibody (orb306039)

Description	Rabbit polyclonal antibody to HIST1H3A
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	H3C1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	15404
Uniprot ID	P68431
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

