

Product Datasheet

TRIM43B antibody (orb306040)

Description	Rabbit polyclonal antibody to TRIM43B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	TRIM43B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52310
Uniprot ID	A6NCK2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

