

Product Datasheet

CYB5R2 antibody (orb306078)

Description	Rabbit polyclonal antibody to CYB5R2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Target	CYB5R2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	31458
Uniprot ID	Q6BCY4
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

