

Product Datasheet

Cbx7 antibody (orb306083)

Description	Rabbit polyclonal antibody to Cbx7
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Target	Cbx7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	18109
Uniprot ID	Q8VDS3
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

