



Product Datasheet

Dlk1 antibody (orb306095)

Description	Rabbit polyclonal antibody to Dlk1
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Target	Dlk1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	41320
Uniprot ID	Q09163
Dilution Range	FACS: 1:25, WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse kidney, m...